

initial data
DICOM files

cohort
folder

patient_1

patient_2

patient_3

1

standardize architecture

Patient_id

Date

Serie number

2

convert dicom to nifti

CT

image.nii.gz

RTSTRUCT

mask_1.nii.gz

mask_2.nii.gz

mask_3.nii.gz

check voxel
intensities

3

tidy files by type

serie
folder

DICOM

NIFTI

4

missing data

- CT serie without RTSTRUCT
- missing nifti present in rtstruct
- tumor mask missing ring mask
- mask baseline missing mask evaluation
- duplicates : same patient id, same date, same location, same ROI num, same ROI type (tum/ring)


Detect problems in you data

5


preprocess nifti

- resample voxel size (mandatory)
- rescale voxel intensity
- apply spatial filter
- discretization of voxel intensity
- delete nifti not resampled

Standardize images and
masks across series and
patients



<https://arxiv.org/html/2402.17521>



<https://forums.ni.com/t5/Machine-Vision/gaussian-filter>





Image intensity discretization. Original data (a) and a generic discretized version (b)

https://www.researchgate.net/figure/mage-intensity-discretization-Original-data-a-and-a-generic-discretized-version-b_fig1_343613289

delete external contours

delete lesions
too small


6 check mask components



decompose tumor masks with several components into a single mask/file per component

delete ring mask with no corresponding ring mask

7 check ring



9 check and correct ring

Create ring when ring incorrect and replace original file

8 create ring

